



ARIA Anesthesia Face Mask *economy without compromise*



- *extraordinarily thin cushion allows for quick and reliable patient seal*
- *non-phthalate patient contact*
- *contains 30% less plastic than a typical anesthesia mask*
- *color coded valve tip for easy size recognition*
- *absolute clarity enables true patient visualization*

U.S. and international patents pending

code	size	dead space (ml)	connector	Qty per box
KM-320	ultra large adult	310	22F	50
KM-321	large adult	250	22F	50
KM-322	medium adult	220	22F	50
KM-323	small adult	180	22F	50
KM-324	child	85	22F	50
KM-325	infant	45	15M	50
KM-326	neonatal	25	15M	50

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY
Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asia and America markets - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong
Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

www.koomedical.co