



www.koomedical.com

## Respiratory Care

### DISPOSABLE MEDICAL FACE MASK

## DISPOSABLE MEDICAL FACE MASK- LARGE

**SYMMETRICAL FOLDING**  
**3D DIMENSIONAL COMFORTABLE**

**BFE≥98%**

Intended use:

The Disposable Medical Face Masks are intended to be worn to protect against the spread or transmission of infectious germs during surgical interventions in operating theatres and other medical facilities. The main aim is to protect the patient against infectious germs. In addition, in certain situations the wearer should be protected against splashes of potentially contaminated liquids and viable particles.

**CLASSIFICATION: TYPE II**

**MODEL: KI-8501**

**STANDARD: EN 14683:2019**

**SPECIFICATION: LARGE**

**SIZE: 17.5cm X 9.5cm**

**QUANTITY: 60PCS**



## DISPOSABLE MEDICAL FACE MASK- SMALL

**SYMMETRICAL FOLDING**  
**3D DIMENSIONAL COMFORT**

**BFE≥98%**

Intended use:

The Disposable Medical Face Masks are intended to be worn to protect against the spread or transmission of infectious germs during surgical interventions in operating theatres and other medical facilities. The main aim is to protect the patient against infectious germs. In addition, in certain situations the wearer should be protected against splashes of potentially contaminated liquids and viable particles.

**CLASSIFICATION: TYPE II**

**MODEL: KI-8502**

**STANDARD: EN 14683:2019**

**SPECIFICATION: SMALL**

**SIZE: 14.5cm X 9.5cm**

**QUANTITY: 60PCS**



### KOO MEDICAL EQUIPMENT

Web: [www.koomedical.com](http://www.koomedical.com)

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: [info@kooeurope.com](mailto:info@kooeurope.com)  
Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23.9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: [sales@koo-asia.com](mailto:sales@koo-asia.com)