




www.koomedical.com

Laryngeal mask

Crystal Airway Mask-Contour Plus

- In transparent PVC- DEHP free.
- Latex-free.
- Soft and non-traumatic laryngeal protector, avoiding lesions and ensuring a good hold and ventilation.
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted.
- Inflation line with valve is integrated in to the wall of the ventilation tube.
- Sterile.
- Single patient use.
- Available in 8 sizes. Each size is color coded.
- Includes gastric channel.



Code	Size	Patient's weight kg	Color	I. D. internal mm	Internal volume of ventilatory pathway (ml)	Items per box
KM-825	5	70~100kg	Yellow	12.8	31	10
KM-835	4	50~70kg	Red	11.8	26	10
KM-845	3	30~50kg	Orange	10.2	16	10