

CPAP Masks

CPAP Zephir



- Safety valve.
- Inflatable cushion with inflation valve.
- It fits perfectly ensuring a secure hold and avoiding air leak.
- Provided with a port for oxygen and monitoring.
- DEHP free.
- PVC is highly transparent and allows good visibility.

- Latex-free.
- Safe and easy to use.
- Provided with two ports: Inhaling (higher) port 22M and exhaling (lower) port (36 mm OD x 33 mm ID).
- Single patient use.
- Available in 3 sizes.

Code	Description	Size	weight g	Dead space ml	Items per box Pcs
KM-301A	mask for CPAP	Large	67.5	245	22
KM-301DA	mask for CPAP	Large deep	74.5	275	22
KM-303A	mask for CPAP	Small	68.5	215	24
KM-303MC	Mask with circuit and manometer for CPAP	Small	-	-	25
KM-301MC	Mask with manometer for CPAP	Large	-	-	25
KM-301DMC	Mask with manometer for CPAP	Large Deep	-	-	25
KM-303C	Mask with circuit for CPAP	Small	-	-	25
KM-301C	Mask with circuit for CPAP	Large	-	-	25
KM-301DC	Mask with circuit for CPAP	Large Deep	-	-	25
KM-301M	Mask with manometer for CPAP	Small	-	-	25
KM-301DM	Mask with manometer for CPAP	Large Deep	-	-	25
KM-303M	Mask with manometer for CPAP	Large	-	-	25

- Each mask is packed individually with a description of the contents.