

Non-Invasive Ventilation Face mask

Bluestar

For lung ventilation devices. Modality CPAP and BILEVEL assisted and/or supervised.



- Inflatable cushion adaptable to the patient's face with a firm and comfortable fit. No skin breakdown.

The kit is provided with:

Transparent PVC mask with inflatable cushion, valve with 6 grips to fix the mask, adapter connector (KM-1019) with ports for O₂ and CO₂, pressure monitoring and a double caps for the introduction of a nasal cannula, DEHP free.

- Pump to inflate and deflate the cushion.
- Headgear's cuff with a rapid removal system.
- Latex Free.
- Possible to link with connector (KM-903 not included) provided with a hole for the CO₂ discharge and a port for the improvement of O₂.
- Single patient use.



The kit is inside a box with a label that indicates the contents.

Code	Description	Size	Weight g	Dead space ml	Dia. of the connector	Items per box Pcs
KPM-001L	face mask NIV	extra large	91.0	290	22F	50
KPM-001M	face mask NIV	large	79.5	240	22F	50
KPM-002M	face mask NIV	medium	70.0	220	22F	50
KPM-003M	face mask NIV	small	65.0	175	22F	50
KPM-004M	face mask NIV	pediatric	46.0	85	22F	50
KM-1019	adapter connector	standard	19.0	17	22M	100