



www.koomedical.com

NIV masks

Non-Invasive Ventilation Face mask

Bluestar PLUS

for pulmonary ventilation's devices, modality CPAP and BILEVEL assisted and supervised.

- Innovative design that guarantees an optimal air-leak control and a high patient comfort level.
- Lightweight.
- Soft replaceable cushion minimizes skin breakdown.
- The kit has all the accessories for positioning.



kit contents:

- Transparent PVC mask with inflatable and replaceable cushion with inflation valve.
- DEHP free.
- Thumbscrew adjustable forehead cushion with two grips.
- Ports for O₂ and CO₂, pressure monitoring and a port for introducing gastric tubing and with the possibility of introducing a bronchoscope.
- Two connectors 22M and 22F.
- Pump for inflating and deflating the cushion.
- Headgear with a rapid removal system.
- Latex free.
- The mask's connector can be connected to expiratory system (KM-903 not in the kit) provided with a hole for the CO₂ discharge and with a port for the enrichment of O₂.
- Every kit is delivered in a box with a label indicating the contents.

Inflatable and replaceable cushions for

Bluestar PLUS

Code	Description	Size	For mask	Items per box
KM-3101	cushion replacement	extra large	KM-310	50
KM-3111	cushion replacement	large	KM-311	50
KM-3121	cushion replacement	medium	KM-312	50
KM-3131	cushion replacement	small	KM-313	50
KM-3141	cushion replacement	pediatric	KM-314	50

Every cushion is packed individually. Mask shells available.



NIV masks



www.koomedical.com

Non-invasive ventilation face mask

Bluestar PLUS



Code	Description	Size	Weight g	Dead Space	Connector	Kits per box Pcs
KM-310B	face mask NIV	extra large	158	205 ml	22F/M	50
KM-311B	face mask NIV	large	150	195 ml	22F/M	50
KM-312B	face mask NIV	medium	140	165 ml	22F/M	50
KM-313B	face mask NIV	small	140	155 ml	22F/M	50
KM-314B	face mask NIV	pediatric	140	95 ml	22F/M	50
KM-310N	Bluestar Plus Kit	Ultra Large	-	205	22F/M	50
KM-311N	Bluestar Plus Kit	Large	-	195	22F/M	50
KM-312N	Bluestar Plus Kit	Medium	-	165	22F/M	50
KM-313N	Bluestar Plus Kit	Small	-	155	22F/M	50
KM-314N	Bluestar Plus Kit	Child	-	95	22F/M	50
KM-310E	kit with AAV-Ultra	Ultra Large	-	205	22F/M	50
KM-311E	kit with AAV	Large	-	195	22F/M	50
KM-312E	kit with AAV	Medium	-	165	22F/M	50
KM-313E	kit with AAV	Small	-	155	22F/M	50
KM-314E	kit with AAV	Child	-	95	22F/M	50
KM-310V	kit with AAV w/o EP	Ultra Large	-	205	22F/M	50
KM-311V	kit with AAV w/o EP	Large	-	195	22F/M	50
KM-312V	kit with AAV w/o EP	Medium	-	165	22F/M	50
KM-313V	kit with AAV w/o EP	Small	-	155	22F/M	50
KM-314V	kit with AAV w/o EP	Child	-	95	22F/M	50

Item KM-310B is provided with mask shell Large and cushion Extra Large.