

For CPAP BILEVEL non-invasive lung ventilation devices



- Transparent PVC mask with inflatable and replaceable cushion with inflation valve.
- DEHP free
- Latex-free
- Designed for comfort and elimination of air leakage.
- New design guarantees even lower skin breakdown.
- Ports with caps for enrichment with O₂ and monitoring.
- Comfortable forehead cushion
- 360 rotating connector
- Single patient use

Delivered in a kit with positioning accessories.

kit contents

• Mask with cushion, inflation pump, headgear.



