



Tri-Gym & Flo-Gym

Lung Exercisers

- In polycarbonate, tube and mouthpiece in polyethylene
- Latex free
- For inspiration and expiration
- · With three measuring cylinders for the control of the flow
- Flow control is done through the three colored balls, which indicate the amount of flow per second 600 cc green, blue 900 cc, 1200 cc red.
- · With corrugated tube and mouthpiece
- · Single patient use



Tri-Gym

Code KM-805

a single packaging, 24 items per box.

Flo-Gym

Can be disassembled, cleaned with mild detergent and reassembled by hand.

Code: KM-802

