Respiratory Care



pipeP Kit

System for lung exercises during inspiration and expiration

- · In polycarbonate, connectors and resistors in ABS
- Latex free DEHP free. Two one-way valves:
 IN= inspiration OD 22mm
 - OUT= expiration OD 22mm ID 15mm
- 8 resistors for PEP or PIP with color code with the following holes' diameters: 1,5 mm black, 2,0 mm white, 2,5 mm yellow, 3,0 mm blue, 3,5 mm green, 4,0 mm orange, 4,5 mm pink, 5,0 mm brown; with connector for manometer, OD 15 mm
- Mask not included in the kit
- Single patient use



Code	Description	Items per box(Pcs)
KM-803	Standard	-
KM-803K	pipeP Connector Kits	90
KM-803HC	pipeP Homecare with Mouthpiece	50
KM-803HCN	pipeP Homecare with silicone mask Child	50
KM-803HCP	pipeP Homecare with silicone mask Child	50
	Large	
KM-803HCA	pipeP Homecare with silicone mask Adult	50
KM-803HCL	pipeP Homecare with silicone mask Extra	50
	Large	