

Evergreen

High concentration oxygen masks

- In transparent PVC
- DEHP free
- Latex-free
- With wide elastic band and adjustable nose clips
- With reservoir and one-way valve
- With 210 cm oxygen supply tubing, kink resistant with a conical standard connector
- Single packaging
- Single patient use



- In transparent PVC.
- DEHP and BPA free
- Latex-free
- With wide elastic band and adjustable nose clips
- VENTURI system with regulator to control the flow of oxygen (24%,26%,28%,31%,35%,40%,50%.)
- With corrugated tube 15 cm
- With 210cm oxygen supply tubing, kink resistant, with a conical standard connector
- Single packaging
- Single patient use

