



www.koomedical.com

# Oxygen therapy

## Evergreen tracheostomy oxygen mask



- In transparent PVC
- DEHP and BPA free
- Latex free
- Fixing elastic band
- Single packaging
- Single patient use

Code	Size	Items per box(Pcs)
KM-908	adult	50
KM-918	pediatric	50

## Evergreen oxygen nasal cannula

- In transparent PVC
- DEHP and BPA free
- Latex free
- With 210cm oxygen supply tubing with a conical standard connector
- Single packaging
- Single patient use

Code	Size	Items per box(Pcs)
with soft and round nasal tip		
KM-218	adult	50
KM-222	pediatric	50
KM-224	infant	50
KM-226	neonatal	50
with soft and flared tip		
KM-228	adult	50
KM-230	pediatric	50

