

## Aerosol FISIOJET kit



- High quality compressor
- Nebulizer Fisiojet KM 1000 with high performance in polycarbonate
- Latex free
- DEHP free
- Elbow connector
- Mouthpiece
- Nose piece
- Connection tube
- In a carry bag

*Every item is in a box*

Code:

FJ001 Junior with Fisiojet infant mask and Fisiojet pediatric mask

FJ002 Senior with Fisiojet adult mask



## FISIOJET nebulizer for aerosol

- High efficiency nebulizer with inspiratory and expiratory valves
- This type of supply optimizes the inhalation therapy, guaranteeing maximum medication absorption during inspiration and minimum medication loss during expiration
- Minimal flow 3.0 l/min maximum flow 6.0 l/m
- Minimum capacity 2 ml maximum capacity 8 ml
- In polycarbonate
- DEHP free
- Latex free

*Provided with the following items:*  
*elbow connector, mouthpiece, nose piece*

- single patient use

Code:

KM-1002 in a box

