

FISIOCHAMBER TRIGON spacer

- **Anti-static Chamber:** the special antistatic treatment of the chamber ensures a complete suspension of the drug, improving the dose available for the patient.
- **Vision Valve:** the movement of the Vision Valve indicates inhalation and permits monitoring of patient compliance. The one-way valve has very little resistance to flow and can be opened with minimal inspiratory effort (suitable for patients with low tidal volume).
- **Universal adapter PROPER FIT:** the soft and silicon material of the adapter ensures a perfect seal with all the inhaler connector available on the market.



- **Spiral vents:** the spiral vent is designed to direct the drug to the one way valve, thereby increasing the quantity of respirable particles
- **Mask COMFORT & SEAL:** the soft material and the ergonomic contour ensures a perfect seal and comfort to the patient in order to allow for a correct inhalation technique. The safety expiratory valve ensures a complete CO₂ exhalation.

Code	Description	Size	Items per box
KM-1030	Chamber	Standard	12
KM-1030B	Chamber, Blue bag	Standard	12
KM-1031	Chamber, Panda Mask	Pediatric Small	12
KM-1031B	Chamber, Panda mask without valve, Blue bag	Pediatric Small	12
KM-1032	Chamber, Panda Mask	Infant	12
KM-1032B	Chamber, Panda mask without valve, Blue bag	Infant	12
KM-1033B	Chamber, Adult Fisio mask, Blue bag	Adult	12

