

Veterinary mask



www.koomedical.com

- Mask for anesthesia, ventilation and aerosol therapy for dogs, cats and small animals.
- Equipped with a black silicone diaphragm and a transparent PVC cone.
- Safe and easy to use.
- Latex-free
- Single packaging



6 sizes available

Code	Size	Weight g	Connector		Items per box Pcs
			External dia.	Internal dia.	
KM-701	6 for big dog	122	27 mm	22 mm	25
KM-702	5 for medium dog	115	27 mm	22 mm	25
KM-703	4 for small dog	47	27 mm	22 mm	25
KM-704	3 for big cat	23.5	27 mm	22 mm	25
KM-705	2 for medium cat	13.5	15 mm	12 mm	25
KM-706	1 for small cat	10	15 mm	12 mm	25

The sizes 3,4,5,6 are delivered with a fixing mask ring

