



www.koomedical.com

# Anaesthesia mask

Anaesthesia Mask

## CRYSTAL



- Inflatable air cushion with inflation valve.
- It fits the patient's face and optimizes the application of anesthetics.
- DEHP free.
- PVC is transparent and allows good visibility.
- Latex-free.
- Safe and easy to use.
- Single patient use.
- Available in 8 sizes with color-coded hook rings.

Code	Description	Size	Weight g	Dead space ml	Colors	Connector	Items per box PCS
KM-200	anaesthesia mask	extra large	59.5	310	transparent	22F	50
KM-201	anaesthesia mask	large	55.5	250	red	22F	50
KM-202	anaesthesia mask	medium	51.0	220	yellow	22F	50
KM-203	anaesthesia mask	small	37.5	180	green	22F	50
KM-204	anaesthesia mask	pediatric	27.0	85	white	22F	50
KM-205	anaesthesia mask	infant	17.5	45		15M	50
KM-206	anaesthesia mask	neonatal	14.5	25		15M	50
KM-207	anaesthesia mask	pre-neonatal	11.0	20		15M	50

- Each mask is packed individually with a description of the contents.
- Appreciated in each country; they guarantee safety and comfort.
- Vanilla scent.

*by request they can be delivered with a strawberry scent with the same code and adding the letter*

*S= KM-200S, KM-201S,  
KM-202S, KM-203S,  
KM-204S, KM-205S,  
KM-206S, KM-207S*

