

Anaesthesia masks CRYSTAL PLUS

- Comfortable inflatable cushion with inflation valve.
- DEHP free.
- PVC is highly transparent and allows good visibility.
- Latex-free
- Safe and easy to use.
- Patented anaesthesia mask designed with lateral grips for a comfortable and confident mask manipulation.
- Single patient use.
- Available in 7 sizes with hook rings.



Code	Description	Size	Weight g	Dead space ml	Connector	Items per box Pcs
KM-279	anaesthesia mask	large	55.5	250	22F	50
KM-278	anaesthesia mask	medium	51.0	220	22F	50
KM-277	anaesthesia mask	small	37.5	180	22F	50
KM-276	anaesthesia mask	pediatric	27.0	85	22F	50
KM-275	anaesthesia mask	infant	17.5	45	15M	50
KM-274	anaesthesia mask	neonatal	14.5	25	15M	50
KM-273	anaesthesia mask	pre-neonatal	11.0	20	15M	50



“CRYSTAL PLUS is considered to be the best mask for anaesthesia. Comfortable and confident manipulation.”

- Each mask is packed individually with a description of the contents.