



www.koomedical.com

Anaesthesia masks

Anaesthesia mask CRYSTAL DIAMOND

- Overmolded cushion technology.
- Comfortable inflatable cushion with inflation valve.
- DEHP free.
- PVC is highly transparent and allows good visibility.
- Latex-free
- Safe and easy to use.
- Mask designed with lateral grips for a comfortable and confident mask manipulation.
- Single patient use.
- Available in 3 sizes with hook rings.



Code	Description	Size	Weight g	Dead space ml	Connector	Items per box Pcs
KM-361	anaesthesia mask	large	48	215	22F	10
KM-362	anaesthesia mask	medium	46	190	22F	10
KM-363	anaesthesia mask	small	38	150	22F	10

- *Each mask is packed individually with a description of the contents.*

