

Crystal Airway Mask



- DEHP free.
- Latex-free.
- Soft and non-traumatic laryngeal protector avoiding lesions and ensuring a good hold and ventilation.
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted.
- Inflation line with valve is integrated in the wall of the ventilation tube.
- Sterile.
- Single patient use.
- Available in 7 sizes. Each size is color coded.

Code	Size	Patient's weight kg	I. D. internal mm	O. D. external mm	Max. volume of the cuff (ml)	Items per box Pcs
KM-877	1	< under 5	5.6	8,6	5	10
KM-867	1.5	5-10	6.4	10	7	10
KM-857	2	10-20	7.8	12	10	10
KM-847	2.5	20-30	9.5	14	15	10
KM-837	3	30-50 small adults and children	10	16	20	10
KM-827	4	50-70 men-women	11.4	18	35	10
KM-817	5	> over 70 adults/men	11.9	20	48	10

