



www.koomedical.com

Reinforced Laryngeal Mask

FLEXCAM

- DEHP free transparent PVC with a metallic reinforcement.
- Latex-free.
- Soft and non-traumatic cuff avoiding lesions and ensuring a good hold and ventilation.
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted.
- Inflation line with valve is integrated in to the wall of the ventilation tube.
- Sterile.
- Single patient use.
- Available in 5 sizes. Each size is color coded.



Code	Size	Patient's weight kg	Items per box Pcs
KM-858	2	10-20	10
KM-848	2.5	20-30	10
KM-838	3	30-50	10
KM-828	4	50-70	10
KM-818	5	> over 70	10

